



Invacare® | **Comet® Ultra**

Built to support with performance  
and reliability



Ideal for those with larger body shapes, the **Comet<sup>®</sup> Ultra** offers a secure, supportive and powerful driving experience without compromising on style or functionality. Packed with all of the safety features and driving characteristics seen throughout the **Comet** family, the **Comet<sup>®</sup> Ultra** will get you to your destination quickly and efficiently, in complete comfort and control.

## Features and options



### Improved suspension

Improved suspension to offer superior comfort to larger individuals.



### Reinforced seat column

Offers additional stability and safety when seated



### Reinforced height and width adjustable armrest

Can be easily adjusted to the user and offers additional support when getting on and off the scooter

## Power, reliability and customisation


















Due to a carefully thought out design process and a robust construction, the safety and stability of the **Comet<sup>®</sup> Ultra** is never a concern; even for users weighing up to 220kg. Developed with the highest quality components and safety features, the **Comet<sup>®</sup> Ultra** makes a trustworthy companion no matter where you go. Furthermore, its flexibility and adaptability means you can customise your scooter to suit your unique style and individual preference.

### Complete comfort and support

The special contoured Ultra Seat has been developed in conjunction with therapists to ensure an optimum seating position and increased comfort. In addition, a pneumatically adjustable lumbar support allows for further comfort when seated if and when required. The reinforced seat base and adjustable armrests offer excellent support and allows the use of the armrest as additional support when getting on and off the scooter.

## Technical data

For more comprehensive pre-sales information about this product, including the product's user manual, please see your local Invacare website.

	Backrest width	Seat depth	Seat to floor height	Backrest height	Driving unit width	Total length	Empty weight	The heaviest part
<b>Comet<sup>®</sup> Ultra</b>	 660 mm	 535 mm	 Min. 440 mm Max. 510 mm	 515 mm	 665 - 685 mm	 1450 mm	 148 kg	 approx. 63 kg
	Max. user weight	Battery capacity	Turning diameter	Max. climbable obstacle height	Ground clearance	Max. safe slope	Max. driving distance	Max. driving speed
<b>Comet<sup>®</sup> Ultra</b>	 220 kg	 2 x 12 V/75 Ah (C20) AGM	 2750 mm	 100 mm	 100 mm	 1 °	 58 km	 10 km/h
	Motor capacity	Colour of shrouds					Upholstery	
<b>Comet<sup>®</sup> Ultra</b>	 10km/h: 650 W / 1800 W <sup>1</sup>	<div> <div>Marble White</div> <div>Sapphire Blue</div> <div>Jade Green</div> <div>Luxor Beige</div> <div>Crystal Copper</div> <div>Black Vinyl</div> </div> <div> <div>Carmine Red</div> <div>Iconic Chrome</div> <div>Diamond Silver</div> <div>Electric Blue</div> <div>Jasper Red</div> </div>						

1. Peak performance

Invacare Australia Pty Ltd:  
ABN 45 074 676 378  
PO Box 5002, 1 Lenton Place, North  
Rocks, NSW 2151, Australia.  
Freephone: 1800 460 460,  
Fax: 02 8839 5353  
Email: [orders@invacare.com.au](mailto:orders@invacare.com.au)  
Web: [www.invacare.com.au](http://www.invacare.com.au)

Invacare New Zealand:  
PO Box 62-124, 4 Westfield  
Place, Mt. Wellington,  
Auckland, New Zealand.  
Freephone 0800 468 222,  
Freefax 0800 807 788  
Email: [sales@invacare.co.nz](mailto:sales@invacare.co.nz)  
Web: [www.invacare.co.nz](http://www.invacare.co.nz)

Invacare Asia Ltd:  
PO Box 5002, 1 Lenton Place,  
North Rocks, NSW 2151,  
Australia.  
Phone: +61 2 8839 5330,  
Fax: +61 2 8839 5331  
Email: [asiasales@invacare.com](mailto:asiasales@invacare.com)  
Web: [www.invacare.com](http://www.invacare.com)

Trusted Solutions,  
Passionate People



© 2017 Invacare International Sàrl. All rights reserved. All information quoted is believed to be correct at time of print. Invacare reserves the right to alter product specifications without prior consultation. Comet Ultra - ANZ - 03/2017