

GO ANYWHERE WITH THESE ALL TERRAIN CHAIRS!!!

SANDCRUISER®

SANDPIPER®



Features:

- All metal components are manufactured from 316L Stainless Steel Metal surface finish electropolished to a mirror surface.
- All tube 1.6 mm wall thickness marine grade Stainless Steel, which means it is tough!!
- Front Casters allow for 360 degree turns
- No special tools required to dismantle for transport, the use of 6 mm drop nose pins & "R" clips make dismantle & reassembly effortless.
- Fully adjustable leg rest with padded hinge protectors.
- Anti tip/brake bar secures the chair in one procedure and also fitted a positive lock / release bar for added safety.
- All assembly pins are secured to the main frame with 40Kg nylon coated stainless steel wire for added strength [note not fitted to rear axle "R" clips]
- Wheelleez™ Low pressure tyres reduces the risk of penetration of sharp objects and produces a very smooth comfortable ride.
- Upholstery is manufactured in "Docril" acrylic canvas which is mould and rot resistant. Docril is also salt water resistant and U/V stable. Seating is also fitted with quick release buckles for easy removal, i.e. cleaning etc.
- A lap belt is fitted to secure the patient to the chair. Additional thoracic belts can be purchased if required.
- Supports friendly...Bodypoint Anterior Trunk/ Aeromesh Chest/Head Supports facilitated plus more.
- Australian Register of Therapeutic Goods (ARTG) No. 191974

Specifications: **SANDPIPER®**

- Seat width 360 mm
- To suit up to 350mm hip width
- Seat height approx. 470 mm
- Length overall 1000 mm including footrest
- Maximum width 840 mm
- Height from ground to top handles 1010 mm
- Front to rear wheels centres 350 mm
- Weight assembled 20 Kg (approx)
- Front: 2 x Wheelleez wheels (24cm x 12 cm)
- Rear: 2 x Wheelleez wheels (42cm x 20cm)



GMS
REHABILITATION

Contact us today for trial and demonstration...

Ph: 1300 734 223 or Email: sales@gmsrehab.com.au

www.gmsrehab.com.au